

Antibody to Melanopsin, AB-N38 Antibody to Melanopsin, affinity-purified, AB-N39

This is an example protocol. Please follow good laboratory technique and safety guidelines.

Working dilutions must be determined for each lot.

Please contact us if you have questions. www.ATSbio.com

IMMUNOSTAINING PROTOCOL

- 1. Remove the corneas, and postfix eyes at 4°C for 24 hours in 4% paraformaldehyde in phosphate-buffered saline (PBS). Remove lenses.
- 2. Cryoprotect eyecups for sectioning at 4°C for 24 hours in 30% sucrose in PBS; embed the eyecups in OCT medium (Sakura Finetek, Torrance, CA), freeze, section (16-20 μ m), and thaw-mount onto gelatin-coated slides.
- 3. Dissect retinas destined for flat-mounting from eyecups immediately after postfixation, stretch onto filter paper, and process in 1.5 ml microfuge tubes.
- 4. Wash tissue (slides and flat-mounts) 3 times (10 min, 4°C) in Tris- buffered saline (TBS, Quality Biological, Gaithersburg, MD) and block for 30 min at 4°C in 1.5% normal goat serum in TBS.
- 5. Incubate tissue for 24 hr at 4°C in a 1:2,500 dilution of antiserum UF006 in a TBS-incubating buffer containing 1% bovine serum albumin, 0.25% carrageenan lambda and 0.003% Triton X-100.
- 6. Wash slides and flat mounts three times in TBS (10 min, 4°C) and incubate for 1 hour at 22°C in Cy3-conjugated anti-rabbit IgG antibody (Jackson ImmunoResearch Laboratories, West Grove, PA) diluted 1:500 in TBS incubating buffer.
- 7. Wash 3 final times in TBS (10 min, 22°C).
- 8. Remove flat-mounts from the filter paper and transfer onto glass microscope slides. Mount flat-mounts and sections in DAPI-containing Vectashield (Vector Laboratories, Burlingame, CA), coverslip, and seal with clear fingernail polish.

Panda S. et al. 2002. Melanopsin (Opn4) requirement for normal light- induced circadian phase shifting Science 298(5601):2213-2216.