



**Anti-Conjugated NO-L-Cysteine  
MOUSE MONOCLONAL  
AB-T125**

**Example of ELISA protocol used to test conjugated NO-L-Cystein:**

1. Coating of conjugated NO-L-Cystein ( $15\mu\text{g/ml}$ ) in maxisorp well plates (Nunc) with a solution of sodium carbonate buffer  $0.05\text{M}$  (pH 9.6), during sixteen hours at  $4^\circ\text{C}$ .
2. Saturation of well plates with of a solution of PBS (pH 7.3) containing  $2.5\text{g/l}$  of BSA (Acros) and  $0.05\%$  Tween 20 (Acros) during one hour at  $37^\circ\text{C}$ .
3. Wash with PBS Tween (two times).
4. Anti-conjugated NO-L-Cystein antibody will be diluted ( $1/1,000$ - $1/5,000$ ) in PBS containing  $2.5\text{g/l}$  BSA and  $10\%$  of glycerol,  $200\mu\text{l}$  by well plate (incubating during 2 hours at  $37^\circ\text{C}$ ).
5. Wash with PBS Tween (three times).
6.  $200\mu\text{l}$  of peroxidase-labeled goat anti-mouse (Sigma) diluted ( $1/10,000$ ) in a solution of PBS containing  $2.5\text{g/l}$  BSA,  $10\%$  of glycerol and  $0.5\%$  of Tween, will be applied by well plate (during one hour at  $37^\circ\text{C}$ ).
7. Well plates will be rinsed with a PBS Tween (three times).
8. And finally the peroxidase will be developed by incubating  $200\mu\text{l}$  by well plate of a citrate  $0.1\text{M}$ /phosphate  $0.2\text{M}$  (pH 5) solution containing  $0.4\%$  of OPD (Sigma) and  $0.03\%$  of hydrogen peroxide (Acros) for ten minutes in the dark, after that, we will stop the reaction by the addition of  $50\mu\text{l}$  of  $2\text{M}$  HCl.
9. The optical density will be measured at  $492\text{nm}$ , to obtain the different values (IC 50).

**Example of Immunohistochemistry protocol**

Perfusion protocol for Adult male Sprague Dawley (weight around 0.5 kg) :

1. The animals can be deeply anesthetized for example with urethane (0.5-1.5g/kg, intraperitoneal).
2. Heparinized, and perfused via the ascending aorta with 100 ml of cold physiologic saline (0.9% NaCl) and with the following fixative solution:
  - a) 300 ml of cold 4% paraformaldehyde and 2% glutaraldehyde in 0.1 M phosphate-buffer (PB), pH 7.2, (twominutes).
  - b) 600 ml of cold 4% paraformaldehyde and 2% glutaraldehyde in 0.1 M phosphate-buffer (PB), pH 7.2, (tenminutes).
  - c) Dissect out the brain and place in a solution of 4% paraformaldehyde in 0.1M PB, pH 7.2, at 4°C for 12-16 hours.
  - d) Before the brains are cut on a freezing microtome, treat the brain in growing concentrations of sucrose (a first solution of 5% of sucrose in PBS until brain sinks), after that, repeat the same process in a solution with a higher level of sucrose (10%), 20%, 25% and finally 30%.

Around 50  $\mu$ m-thick serial sections should be obtained, kept at 4°C in PBS (0.1 M, pH 7.2) and processed for immunostaining.

**Example of immunohistochemical protocol :**

1. In order to avoid possible interference with endogenous peroxidase, free-floating sections are treated with distilled water containing NH<sub>3</sub> (20%), H<sub>2</sub>O<sub>2</sub> (30%) and NaOH (1%) for 20 min (another method is using a solution with 33% of H<sub>2</sub>O<sub>2</sub> and 66% of methanol).
2. Wash the sections for 20 min in 0.15 M phosphate-buffered saline (PBS) (pH 7.2)
3. Pre-incubate for 30 min in PBS containing 10% normal horse serum and 0.3% of Triton X-100 (mixed solution).
4. Incubate at room temperature (1h 30min) and overnight at 4°C in the same mixed solution containing anti-conjugated NO- L-Cysteine antibody (diluted 1/1,000 to 1/5,000; as recommended dilution).
5. Wash the sections in PBS (30 min).
6. Incubate for 60 min at room temperature with biotinylated anti-mouse immunoglobulin (Vector) diluted 1/200 in PBS.
7. Wash the sections in PBS (30 min).
8. Incubate sections for 1 h with a 1/100 diluted avidin-biotin-peroxidase complex (Vectastain).
9. Wash the sections in PBS (30 min).
10. Wash with Tris-HCl buffer, pH 7.6 (10 min).
11. Develop tissue-bound peroxidase with H<sub>2</sub>O<sub>2</sub> using 3, 3' diaminobenzidine as chromogen.
12. Rinse sections with PBS and coverslip with PBS/Glycerol (1/1).



**Anti-Conjugated NO-L-Cysteine  
MOUSE MONOCLONAL  
AB-T125**

**Example of Western blot protocol**

Membrane Blocking, Antibody Incubations and Detection

1. Saturate the blot membrane with TBS + 5% Blocker for 1 hour at 37°C while mixing
2. Wash the membrane twice for 5 minutes in TBS Tween at 37°C
3. Incubate the membrane with anti-NO-L-Cysteine antibody diluted 1:1000 in TBS 0.5% Blocker for 2 hours at 37°C
4. Wash the membrane three times for 5 minutes in TBS Tween at 37°C
5. Incubate with a biotinylated secondary antibody diluted (1/1,000-1/2,000) in TBS 0.5% Blocker for 2 hours at 37°C
6. Wash the membrane three times for 5 minutes in TBS Tween at 37°C
7. Incubate with Streptavidin-HRP 1µg/ml in TBS 0.5% Blocker for 2 hours at room temperature
8. Wash the membrane three times for 5 minutes in TBS at 37°C
9. Incubate in TBS (200ml) + (50mg DAB in 25ml methanol) + (150mg 4-chloro-1-naphtol in 25ml methanol) + 50µl H<sub>2</sub>O<sub>2</sub> 30% for a maximum of 30 minutes in the dark
10. Stop the reaction by addition of distilled water

Blocker = skim milk (Biorad 170-6404)

TBS = 20mM Tris base, 0.5M NaCl, pH 7.5

TBS Tween = TBS + 0.05% Tween 20