



**Anti-Conjugated Ascorbic Acid  
RAT POLYCLONAL  
AB-T151**

**Example of ELISA protocol used to test conjugated ascorbic acid:**

1. Coating of conjugated ascorbic acid ( $15\mu\text{g/ml}$ ) in maxisorp well plates (Nunc) with a solution of sodium carbonate buffer 0.05M (pH 9.6), during sixteen hours at  $4^{\circ}\text{C}$ .
2. Saturation of well plates with of a solution of PBS (pH 7.3) containing 1g/l of BSA (Acros) and 0.05% Tween 20 (Acros) during one hour at  $37^{\circ}\text{C}$ .
3. Wash with PBS Tween (three times).
4. Preabsorbed ascorbic acid antiserum will be diluted (1/2,000-1/5,000) in PBS tween containing 1g/l BSA and 5% of glycerol,  $200\mu\text{l}$  by well plate (incubating during 2 hours at  $37^{\circ}\text{C}$ ).
5. Wash with PBS Tween (three times).
6.  $200\mu\text{l}$  of peroxidase-labeled goat anti-rat (Jackson) diluted (1/5,000) in a solution of PBS containing 1g/l of BSA and 0.5% of Tween, will be applied by well plate (during one hour at  $37^{\circ}\text{C}$ ).
7. Well plates will be rinsed with a PBS Tween.
8. And finally the peroxidase will be developed by incubating  $200\mu\text{l}$  by well plate of a citrate 0.1M/phosphate 0.2M (pH 5) solution containing 0.4% of OPD (Sigma) and 0.03% of hydrogen peroxide (Acros) for ten minutes in the dark, after that, we will stop the reaction by the addition of  $50\mu\text{l}$  of 2M HCl.
9. The optical density will be measured at 492nm, to obtain the different values (IC 50).