

Targeting Talk: Agonist Effects

by Dr. Douglas Lappi

Q: Do conjugated toxins (Dermorphin-saporin in particular) exhibit agonist effects? I've generated behavioral and tissue time course effects but have not established agonist effects for this conjugated toxin.

A: The peptide ligand toxins should exhibit agonist effects. They are constructed purposely to retain complete agonist activity, including for us the most important: internalization. So, for instance, SP-SAP (Cat. #IT-07) causes receptor internalization similar to SP, as reported in Mantyh *et al.* (Inhibition of hyperalgesia by ablation of lamina I spinal neurons expressing the substance P receptor. *Science* 278:275-279, 1997).

As to dermorphin-SAP (Cat. #IT-12) specifically, it has agonist activity very much like dermorphin. This is reported in Porreca *et al.*

(Inhibition of neuropathic pain by selective ablation of brainstem medullary cells expressing the mu-opioid receptor. *J Neurosci* 21:5281-5288, 2001) in which it's stated:

The bilateral microinjection of 3 pmol of dermorphin or of dermorphin-saporin directly into the RVM produced a robust antinociceptive effect in the 52°C hot-water tail-flick test. The peak antinociceptive effect of dermorphin, $78 \pm 13.2\%$ MPE, was not significantly different from that of the dermorphin-saporin conjugate, which was $59 \pm 4.7\%$ MPE ($p > 0.5$, Student's *t* test).

Usually the amount needed to give a response is lower than the amount needed to kill a cell. Depending on what your system is; it may be a peculiarity of that system, but I would be a little concerned about not seeing an agonist effect. On the other hand, if you have demonstration of specific toxicity, it may not be all that crucial.

SP-SAP Research Replacement

SP-SAP (formerly Cat. #IT-07) is being replaced by SSP-SAP (Cat. #IT-11). Scientific advisors have given counsel for this replacement because SSP-SAP is a superior lesioning agent in many situations, due to its targeting vehicle, a protease-resistant form of substance P. An excellent example is the paper by Martin and Sloviter, *J Comp Neurol* 436:127-152 (2001), in which after ineffective intraparenchymal injection of SP-SAP in the hippocampus, SSP-SAP was used with tremendous efficacy. In almost all applications we expect SSP-SAP will be used at a lower dose than SP-SAP.

We will not compromise the work of those that are midstream in experiments with SP-SAP. Orders from those who wish to continue with SP-SAP will be honored. This will be a gradual phase-out since SP-SAP is being developed as a chronic pain therapeutic (see update on Page 2). Please rest assured that your work is the backbone of our company and all questions or comments will be carefully considered.

**Questions about an ATS product
or Molecular Surgery technique?
Write to us at ats@ATSbio.com**

Flow Cytometry Services

Cell Sorting \$75/hr

Cell Analysis

**Cell Cycle, Apoptosis,
Antibody specificity**

**1-10 samples @ \$15 ea
>10 samples ask for quote**

Lab Services ask for quote
**Antibody titration, cell culture,
Staining, Blood processing**

(858) 642-1988

Cytometry Research, LLC

11175-A Flintkote Ave, SD 92121

www.CytometryRes.com

A subsidiary of Advanced Targeting Systems