

Targeting Tools: Featured Products


mu p75-SAP

In 2004, ATS re-designed the anti-murine p75-SAP targeted toxin (mu p75-SAP, Cat. #IT-16) and produced a conjugate that is much more potent in our *in vitro* cell cytotoxicity assays. Previously, we used a rat monoclonal antibody. This antibody had been outperformed by our rabbit polyclonal (Cat. #AB-N01), in several assays, especially flow cytometry analysis of murine p75-expressing cells. This is an important indicator of being able to bind to the cell surface, which is fundamental for a targeted toxin.

To create this toxin, we affinity-purified the rabbit polyclonal (Cat. #AB-N01AP) with the immunogen bound to a solid support, and conjugated the affinity-purified antibody to saporin. As can be seen in the cytotoxicity assay on the right, the new mu p75-SAP is orders of magnitude more potent than the previous conjugate. The new and more active version of mu p75-SAP has an ED₅₀ in the picomolar range compared to an ED₅₀ in the nanomolar range for the previous product. We believe that the greater potency will translate to smaller amounts used for elimination of p75-positive neurons in the mouse brain, and that this will result in a greater index of efficacy and lesser non-specific cytotoxicity. (see cover article for *in vivo* results).

The mu p75-SAP kit includes, in addition to the immunotoxin, equal aliquots of saporin (Cat. #PR-01), the affinity-purified rabbit polyclonal antibody (AB-N01AP), and the control immunotoxin, Rabbit-IgG-SAP (Cat. #IT-35).

Also available are fluorescent conjugates of AB-N01AP: Cy3-labeled Anti-murine NGFr (Cat. #FL-05), and Cy5-labeled Anti-murine NGFr (Cat. #FL-06).



NG3 cells are plated at 1000 cells/well and incubated overnight. Saporin, mu p75-SAP (conjugate of the affinity-purified rabbit polyclonal to mouse NGFr and saporin), and AB-N02-SAP (previous rat monoclonal version of mu p75-SAP) are added in 10- μ l volumes and the plates are incubated 72 hours. PMS/MTS developing reagent is added and the plates are incubated 1-2 hours, then read at 490 nm.

Flow Cytometry Tools

Advanced Targeting Systems is excited to offer two new tools for flow cytometry.

WBlyse™ (Cat. #FL-08) is a gentle erythrocyte lysing reagent matched with a leukocyte preservative. This kit can be used to enumerate lymphocyte subsets, detect a large variety of antigens on lymphocytes, and to identify other leukocyte subsets, including CD34 stem cells, and granulocytes. WBlyse™ works on many types of specimens and is active against all erythrocytes and is available in the 100 or 500 test size.

The CFSE (carboxyfluorescein-succinimidyl ester) Compensation Kit (Cat. #FL-13) provides a reproducible and convenient source of CFSE particles for evaluating and correcting for the spectral overlap between CFSE and other dyes. Fluorescence compensation is the process wherein a portion of the primary dye signal is removed from all non-primary dye channels. Compensation is generally required whenever more than a single dye (color) is used in flow cytometric analysis. This 1-color CFSE Compensation Kit is sufficient for 25 tests and contains negative/unstained microspheres and CFSE microspheres. Cells, once labeled with CFSE, can be stimulated *in vitro* and cell proliferation measured by changes in the staining intensity of the cells.

