

Saporin Uses in Immunology

One of the more common methods of using Streptavidin-ZAP is to couple the complex with biotinylated antibodies. However, there are also many instances of biotinylated ligands, peptides, and even biotinylated major histocompatibility complex (MHC) tetramers in immunology fields.

B-Cell Targeting. An example of utilizing Streptavidin-ZAP to study ligands that target CD22 were shown by Collins et al. [1]. CD22 is a potential target for immunotherapy of B-cell lymphomas. The authors examined the equilibrium between CD22 and the cis and trans forms of its ligands using high affinity sialoside probes. They also demonstrated that a biotinylated probe specific for CD22, when combined with Streptavidin-ZAP, can eliminate several different lymphoma cell lines.




Figure 5. Selective cytotoxicity by enzymatic activation of CCP1-SA-ZAP conjugates. Schematic representation of streptavidin-ZAP bound to the CXP1 peptides (A) and the expected toxicity of the different SA-ZAP conjugates to ACPA expressing B cells (B). Percentage of living ACPA-expressing B cells (C) and TT-specific B cells (D) after 4 days of treatment with antigen-toxin conjugates [2].

Rheumatoid arthritis is a chronic disease that is accompanied with anti-citrullinated protein antibodies (ACPA) produced by autoreactive B cells. A study in 2018 used a synthesized cyclic citrullinated peptide (CCP) antigen suitable for B-cell receptor binding and demonstrated that binding by ACPA was impaired upon manipulation of the residue [2]. The data were generated using biotinylated CCP mixed with Streptavidin-ZAP in cell viability assays. The results marked an important step towards antigen-selective B-cell targeting in general, and more specifically in rheumatoid arthritis.

T-Cell Targeting. An example of using Streptavidin-ZAP to deplete specific T cells was published in 2010; Akiyosi et al. used the dendritic cell-associated heparan sulfate proteoglycan-dependent integrin ligand (DC-HIL) as the targeting agent. DC-HIL is exclusively associated with syndecan-4 (SD-4) which is expressed on some activated T cells [3]. A similar study was done with Sézary syndrome cells that overexpress syndecan-4 [4].

Hess et al. investigated whether pathogenic T cells could be depleted via Streptavidin-ZAP coupled to MHC class I tetramers to kill antigen-specific CD8+ T cells [5]. Their work showed the therapeutic potential for using cytotoxic tetramers to eliminate specific T cells. This same strategy was employed *in vivo* to delay diabetes in non-obese diabetic mice [6]. The Hess group also used biotinylated peptide-MHC class I tetramers with Streptavidin-ZAP to selectively deplete a population of alloreactive T cells in mice to determine that toxic tetramer administration prior to immunization increased survival of cognate peptide-pulsed cells in an *in vivo* cytotoxic T lymphocyte assay and reduced the frequency of corresponding T cells [7]. More research towards T-cell depletion utilizing Streptavidin-ZAP and biotinylated MHC tetramers came from Sims et al. where they found that following a significant transient depletion of cells, the population rebounded and reached a higher percentage of total CD8+ T cells than before the depletion. This research provides helps further understanding of the ‘flexibility and turnover’ of these cells [8].




Figure 1. Frequency and function of MCMV-specific CD8+ T cells and depletion of M38-specific CD8+ T cells. C57BL/6 mice were infected *i.v.* with 1×10^6 pfu MCMV. **D)** Representative flow cytometry plots of M38-specific CD8+ T cells and M45-specific CD8+ T cells 50 days after MCMV infection, 1 and 6 days after either M38-tetramer-saporin (red box) or M38-tetramer-PE injection (orange box). Data are shown as mean \pm SEM and are from one representative plot out of 6/group) [8].




Figure 7. Sequential treatment with mAb 440c and saporin depletes IPCs in vitro. (A) Expression of 440c on bone marrow cells cultured in FLT3-L for 10 days. CD11c⁻ cells were excluded from the gate. 440c is highly expressed on B220⁺ cells, which represent fully developed IPCs. (B-C) Depletion of CD11c⁺/B220⁺ cells from FLT3-L-derived bone marrow IPCs by treatment with 440c-saporin. Cells were stained on ice with biotinylated-440c (B) or a biotinylated control antibody (C), followed by avidin-saporin, cultured for an additional 36 hours, and analyzed for residual presence of CD11c⁺/B220⁺ cells. The percentage of gated cells is indicated. Results are representative of 3 separate experiments [10].

Other Applications. Outside of T cells and B cells, Streptavidin-ZAP has also been used to deplete dendritic cells (DC). Alonso et al. depleted inflammatory DCs with biotinylated anti-CD209 via intravenous injection in a mouse animal model of induced inflammatory DC formation [9]. The authors suggest that the depletion of inflammatory DCs could be useful in understanding inflammatory diseases such as psoriasis. Depletion of natural interferon-producing cells (IPC) was demonstrated with an IPC-specific biotinylated antibody in vitro [10].




Figure 3. CCL28 promotes tumour growth through attracting CCR10⁺ Treg cells. **h)** CCR10 depletion eliminates most CD4⁺CD25⁺FOXP3⁺ cells but not CD8⁺ T cells, as determined by flow cytometric analysis. CCR3 depletion eliminates both populations (numbers inside plots refer to the boxed areas: left column, % Treg cells; right column, % CD31⁺CD41⁺ cells (left) and % CD31⁺CD81⁺ cells (right)). or M38-tetramer-PE injection (orange box). Data are shown as mean ± SEM and are from one representative plot out of 6/group) [8].

Facciabene et al. made immunotoxins with Streptavidin-ZAP to biotinylated antibodies for CCR10 and CCR3 (Anti-CCR10-Saporin and Anti-CCR3-Saporin). They investigated whether a direct link between tumor hypoxia and tolerance occurs through the recruitment of regulatory cells [11]. Their findings showed that peripheral immune tolerance and angiogenesis programs are closely connected and cooperate to sustain tumors.

References

1. Collins, B.E.; Blixt, O.; Han, S.; Duong, B.; Li, H.; Nathan, J.K.; Bovin, N.; Paulson, J.C. *J Immunol* **2006**, *177*, 2994-3003, doi:10.4049/jimmunol.177.5.2994.
2. Lelieveldt, L.; Kristyanto, H.; Pruijn, G.J.M.; Scherer, H.U.; Toes, R.E.M.; Bongers, K.M. *Mol Pharm* **2018**, *15*, 5565-5573, doi:10.1021/acs.molpharmaceut.8b00741.
3. Akiyoshi, H.; Chung, J.S.; Tomihari, M.; Cruz, P.D., Jr.; Ariizumi, K. *J Immunol* **2010**, *184*, 3554-3561.
[Check out the Targeting Trends article on our website.](#)
4. Chung, J.S.; Shiue, L.H.; Duvic, M.; Pandya, A.; Cruz, P.D., Jr.; Ariizumi, K. *Blood* **2011**, *117*, 3382-3390.
5. Hess, P.R.; Barnes C Fau - Woolard, M.D.; Woolard Md Fau - Johnson, M.D.L.; Johnson Md Fau - Cullen, J.M.; Cullen Jm Fau - Collins, E.J.; Collins Ej Fau - Frelinger, J.A.; Frelinger, J.A. *Blood* **2007**, *109*, 3300-3307, doi:10.1182/blood-2006-06-028001.
[Check out the Targeting Trends article on our website.](#)
6. Vincent, B.G.; Young, E.F.; Buntzman, A.S.; Stevens, R.; Kepler, T.B.; Tisch, R.M.; Frelinger, J.A.; Hess, P.R. *J Immunol* **2010**, *184*, 4196-4204, doi:10.4049/jimmunol.0903931.
7. Hess, S.M.; Young, E.F.; Miller, K.R.; Vincent, B.G.; Buntzman, A.S.; Collins, E.J.; Frelinger, J.A.; Hess, P.R. *Transl Immunol* **2013**, *29*, 138-145, doi:10.1016/j.trim.2013.10.005.
8. Sims, S.; Bolinger, B.; Klenerman, P. *Eur J Immunol* **2015**, *45*, 113-118, doi:10.1002/eji.201445016.
9. Alonso, M.N.; Gregorio, J.G.; Davidson, M.G.; Gonzalez, J.C.; Engleman, E.G.. *Immunol Res* **2014**, *58*, 374-377, doi:10.1007/s12026-014-8511-6.
10. Blasius, A.; Vermi, W.; Krug, A.; Facchetti, F.; Cella, M.; Colonna, M. *Blood* **2004**, *103*, 4201-4206, doi:10.1182/blood-2003-09-3108.
11. Facciabene, A.; Peng, X.; Hagemann, I.S.; Balint, K.; Barchetti, A.; Wang, L.P.; Gimotty, P.A.; Gilks, C.B.; Lal, P.; Zhang, L.; et al. Tumour hypoxia promotes tolerance and angiogenesis via CCL28 and T(reg) cells. *Nature* **2011**, *475*, 226-230, doi:10.1038/nature10169.