



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

**HRP-labeled Antibody to HA
MOUSE MONOCLONAL**

Catalog Number: AB-466
Quantity: 500 microliters
Format: 1×PBS and 50% glycerol
Host: Mouse
Isotype: IgG2b
Clone: PHAHRPSHG
Immunogen: Synthetic peptide (YPYDVPDYA) coupled to KLH

Specificity and Preparation:

Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (YPYDVPDYA) conjugated to KLH. The antibody is conjugated to horseradish peroxidase (HRP).

Usage and Storage:

Reported to be effective for immunoblotting (western blot, 1:500).

Antibody is stable for two years at -20°C, one month at 4°C. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols