



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

Antibody to HA
RABBIT POLYCLONAL

Catalog Number: AB-469
Quantity: 100 microliters, 200 microliters, 500 microliters
Format: Serum with 0.02% sodium azide
Host: Rabbit
Clone: PPHASHG
Immunogen: Synthetic peptide (YPYDVPDYA) coupled to KLH

Specificity and Preparation:

Polyclonal antibodies are produced by immunizing rabbits with synthetic peptide (YPYDVPDYA) conjugated to KLH.

Usage and Storage:

Reported to be effective for immunoblotting (western blot, 1:1,000) and immunoprecipitation. Antibody is stable for two years at -20°C, one month at 4°C. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols