



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

Dengue Virus NS1 Recombinant
VIRAL ANTIGEN

Catalog Number: PRP-404
Quantity: 100 micrograms, 500 micrograms, 1 milligram
Format: 100mM Tris HCl pH-8, 8M Urea & 200mM NaCl
Host: *E. coli*

Background:

Dengue virus is caused by one of four closely related virus serotypes of the genus *Flavivirus*, family Flaviviridae. Each serotype is sufficiently different from the other such that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and in mice Morpholino antisense oligos have shown specific activity against dengue virus.

Specificity and Preparation:

The *E. coli*-derived recombinant protein contains the NS1 dengue virus full length type-2 immunodominant regions. The dengue protein is fused to 6His tag. It is purified by proprietary chromatographic technique. The protein is >90% pure as determined by PAGE followed by coomassie staining.

Usage and Storage:

The end user must determine the optimum working titer for each particular application. Although stable at 4°C for 1 week, the material should be stored below -18°C. Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols